**PROTOCOL FOR USING THE HEMOCUE™**

**ENGLISH/SPANISH VERSIONS**

## **MOTHERCARE USING THE HEMOCUE™<sup>1</sup>**

Many community-based health programs use the HemoCue™ to assess anemia prevalence, based on low hemoglobin. The purpose of this paper is to provide information on the accuracy of the HemoCue™ and how to maximize its accuracy in the field setting. Information from the U.S. Centers for Disease Control and Prevention (CDC) and HemoCue Inc. was used to prepare this paper. A step-by-step procedure for using the HemoCue™ in programs is attached.

Staff at CDC feels that the HemoCue™ is a robust instrument that can give accurate readings in a field setting. Its technology is much more sophisticated than the older hemoglobinometer previously relied on. However, there are some subtle techniques in using this instrument that, if not followed, may cause wide variation in results. The greatest variation arises when capillary blood is not sampled properly. If an adequate puncture (finger-stick or toe-stick) is not made to allow for spontaneous blood flow, the mix of red blood cells in the sample will not be representative (red blood cells concentration may be too low or too high). The hemoglobin (Hb) determination will not accurately represent the patient's blood hemoglobin level and could lead to false-positive or false-negative determination of anemia. The issue of a representative blood sample is not a problem when venous blood is obtained.

Wide variations in Hb results can have significant impact on determination of anemia prevalence. Based on data from a national survey in the United States (NHANES II), the expected standard deviation of the Hb distribution for non-iron deficient one to five year old children is 0.84 using venous blood. A wider standard deviation will increase anemia prevalence even if mean hemoglobin does not change. For example, if the standard deviation is increased to 1.2, the anemia prevalence changes from 9.5% to 17.9%.

---

<sup>1</sup>This paper was prepared by Rae Galloway, MotherCare Nutrition Advisor (1994-1998) with information and comments from Caryn Bern and Ibrahim Parvanta at the Centers for Disease Control and Prevention, Atlanta, Georgia.

The standard deviation will be higher for capillary sampling. It also will be higher for pregnant women because of hemodilution, particularly if they are sampled throughout pregnancy. While it is not known what a reasonable variation should be for pregnant women, specialists<sup>2</sup> at CDC suggest that a standard deviation of 1.2 is reasonable for pregnant women if they are sampled from mid-second trimester to parturition.

From work which CDC has conducted in the Supplemental Nutrition Program for Women Infant and Children (WIC) clinics in the United States, training can improve health workers' techniques for collecting capillary blood using the HemoCue™. Staff members at CDC feel that in a research setting where one-time anemia prevalence data are collected, training is useful for decreasing large variation in hemoglobin results. In a program setting where the HemoCue™ is used for screening purposes, training will need to be reinforced periodically and health workers supervised to maintain high quality capillary blood sampling techniques and data reporting.

The following are the most common problems in capillary sampling and the use of the HemoCue™:

**1) Keep the instrument clean. This applies especially to the cuvette holder.**

A swab dabbed with alcohol can be used to clean away any dirt or dried blood. This should be done at least once a day or when there is a visible build-up of dirt or blood.

**2) Ensure instrument accuracy.**

Check the accuracy of the instrument daily, or when performance is questioned, using the control cuvette that comes with each HemoCue™ instrument. Keep a daily log of accuracy readings. If the accuracy readings are outside the range of the control cuvette, and the HemoCue™ is clean, then the instrument needs to be calibrated by a trained technician.

**3) Make sure the finger stick is adequate.**

Wide variations can occur in Hb measurements from the same person if the finger stick is inadequate (i.e., the finger stick is not deep enough to allow adequate flow of blood and a representative concentration of red blood cells). In most cases if the finger stick is inadequate, hemoglobin (Hb) values will be underestimated and the

---

<sup>2</sup>Ray Yip and Ibrahim Parvanta, personal communication .

prevalence of anemia will be overestimated. The type of lance used can mitigate this error. It should be noted that in all MotherCare programs a retractable lance should be used to ensure it is not re-used.

#### 4) **Avoid poor technique.**

- a. Milking the finger<sup>3</sup> (usually related to an inadequate finger stick) to obtain proper blood flow which will underestimate Hb and overestimate anemia prevalence.
- b. Lack of uniformity as to which drop of blood is taken. Wipe away the first three drops of blood using dry wipe and measure the fourth or subsequent drops [HemoCue™ recommends this]. This helps allow for the collection of blood with a representative mix of red blood cells.
- c. Mixing the blood with alcohol. The finger should be totally dried before the finger prick is performed. Use alcohol before the prick is made and then wipe away each drop of blood with a dry wipe to avoid any mixing of the blood with alcohol. Wiping away the first few drops of blood also will minimize the mixing of sweat with blood in hot, humid climates. This error will underestimate Hb and overestimate anemia prevalence.

#### 5) **Adequately fill the microcuvette.**

The microcuvette needs to be filled with an entire drop of blood in one attempt. Again this depends on the flow of blood; if it is not adequate, the microcuvette will not fill adequately. Topping off the microcuvette with more than one attempt will overestimate the Hb reading and underestimate the anemia prevalence.

Any signs of air-bubbles means that the microcuvette has not been filled adequately and should be discarded and a new microcuvette used. The presence of bubbles will underestimate the Hb reading and overestimate anemia prevalence.

#### 6) **Take the reading within 10 minutes after collecting the blood sample.**

Delays in reading will result in erroneous Hb values.

---

<sup>3</sup>When the finger is squeezed to increase blood flow.

## Step-by-Step Use of the HemoCue™

Each project should assign one person who, over a period of a few days, will repeatedly utilize the HemoCue™ to decrease the variability in his/her technique. This process needs to be documented and related in the methodology section of each baseline survey report. The trainer should ensure that trainees obtain the proper technique to give them a standard deviation of 1.2. It is suggested that trainees practice on each other to acquire the proper technique.

### Daily

At the beginning of each day, check instrument accuracy using the control cuvette. If readings are in question, clean the cuvette holder and check the accuracy again. If readings continue to be outside the correct range, the instrument should be re-calibrated by a trained technician. (In the MotherCare project sites, this means that the instrument should be sent back to HemoCue™ for adjustment since we will not have the capacity in the field to re-calibrate using a blood standard and making an internal adjustment.)

Check microcuvette containers for the expiration date that is good if the container is unopened (sealed). Opened microcuvette containers are good for 3 months after opening. Upon opening, the container should be clearly dated with the open date. Keep the container lid closed when not use to avoid unnecessary exposure of the cuvettes to air, especially in humid conditions. Heat and moisture will denature the chemicals in the cuvette that can lead to inaccurate Hb measurements.

### For each use:

1. Check the instrument for dirt or dried blood -- particularly the microcuvette holder -- and use a cotton swab dipped in alcohol to clean the affected area. Allow this to air dry completely.
2. Seat the patient comfortably. Ensure that the patient's hand is warm so that blood circulates freely before sampling. The patient's fingers should be straight but relaxed to avoid the stasis effect that occurs when the fingers are bent.
3. Use only the middle finger or ring finger for sampling. Remove rings from the finger. Clean the puncture site with a disinfectant and allow it to dry.
4. Using a rolling movement of your thumb, lightly press the finger from the top knuckle towards the tip. This stimulates the flow of blood towards the sampling area.

5. When the thumb has reached the fingertip, maintain gentle pressure, and prick the side of the tip with a sharp, quick motion. This causes the least pain and the best flow of blood. Use a retractable lance and dispose of it in the appropriate receptacle.
6. Using a dry absorbent pad, wipe away the first three drops of blood to stimulate spontaneous blood flow. If necessary, press gently again until another drop of blood appears. Avoid "milking."
7. Fill the microcuvette by touching its pointed tip to the middle of the blood drop. This helps avoid trapping air in the microcuvette. Ensure that the drop of blood is big enough to fill the microcuvette completely in one attempt (the microcuvette fills automatically by capillary action). Never top up the microcuvette after the first filling.
8. If a second sample is taken from the same site, wipe away the remains of first sample, apply gentle pressure to form another adequate blood drop and collect the sample as described above.
9. Wipe off any excess blood on the sides of the microcuvette. Ensure that no blood is sucked out of the microcuvette.
10. Place the microcuvette in its holder and gently push the holder into the photometer. The microcuvette should be read within ten minutes (at the latest) after being filled. It takes 14-45 seconds for a reading. Reseal the microcuvette container immediately after use.

### **Cleaning the Hemocue™**

Cleaning the photometer is probably the most effective method for preventing and solving problems. Following are instructions to most efficiently and effectively clean the instrument. (Materials needed include cotton swabs and water).

#### **Cleaning the microcuvette holder:**

Remove the black microcuvette holder and clean it with a gentle cleanser and hot water or alcohol. Rinse it completely with water and dry it well.

#### **Cleaning the Optronic unit of the photometer**

- a. Wet the cotton swab with water and squeeze out excess water.

- b. Having already taken out the cuvette holder, insert the swab about 5 centimeters in the cavity that holds the cuvette holder. By applying pressure to the roof of the cavity, you will feel another small cavity. This is the sensor window.
- c. Clean an area of about one inch surrounding the window. Repeat this process until the swab no longer picks up any RESIDUAL blood.
- d. Using another swab, go over the area you just cleaned until it is completely dry.
- e. Ensure the proper operation of the photometer by using the control cuvette.

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